

PROTEOR QUATTRO[®] PATIENT TRIAL / DELIVERY CHECKLIST

GENERAL INFORMATION

Date: _____

Facility: _____

CP: _____

Patient ID: _____

Pyramid or Threaded Top: P T

Foot: _____ Side: R L Size: _____ Cat: _____ Toe: Sandal Reg Footshell Color: L M D B

FACILITY CHECKLIST

- Access to Ramps or Slopes
- Access to Stairs
- Access to Uneven Terrain

FOR DELIVERY START AT STEP 3

STEP 1: SET EXPECTATIONS (PATIENT AND CP)

- The QUATTRO is a revolutionary design and will not walk like your current knee
- Discuss how the valve control technology is constantly making changes to optimize gait on every step or situation
- Consider taking a video of each knee to provide a comparison (Consent form)

STEP 2: ACTIVITIES ON CURRENT KNEE

Prior to the patient trialing the QUATTRO, have the patient perform these tasks on their current knee:

- Walking Backwards
- Side steps
- Slow, Medium, Fast walking
- Transition from slow to fast, fast to slow
- Stairs and Ramps step over step (if they currently do)
- Sitting and sit to stand
- Uneven Terrain (if available)

STEP 3: INSTALL AND ALIGN THE QUATTRO

- Ensure proper alignment prior to programming the knee
- Bench Alignment (no Stance Flexion: 0-5mm anterior; with Stance Flexion 0-5mm posterior to knee center)
- Static Alignment- patient standing quietly (no Stance Flexion: 0-10mm anterior; with Stance Flexion 0-10mm posterior to knee center)
- Dynamic Alignment- patient walking

STEP 4: COMPLETE PROGRAMMING IN GAITLAB

- Calibration (consider walking outside the parallel bars with walking aid for Step 3)
- Resistance Settings
- Demonstrations

PROTEOR QUATTRO[®] PATIENT TRIAL / DELIVERY CHECKLIST

STEP 5: DEMONSTRATIONS HIGHLIGHTING FEATURES

- Walking Backwards
- Side steps
- Slow, Medium, Fast walking
- Transition from slow to fast, fast to slow
- Stairs and Ramps step over step (if they currently do)
- Sitting and sit to stand
- Uneven Terrain (if available)

	KEYPAD	GAITLAB	FI APP
Stance Lock (on/off if applicable)	—		
Obstacle Assist	—		—
Flexion Lock			
Sitting Mode	—		
Select Activities	—		
Setting Alternate Activities	—		
Switch Activities (Walking and Alternate)		—	—
Add New Activities	—		
Adjust Activities	—		
Cycling	—		
Golfing	—		
Sprinting	—		
Elliptical	—		
Rehabilitation	—		
Remote Firmware Update	—		
Cadence and Usage Reports	—		
Outcome Measures	—		

FINAL SETTINGS

Sitting: _____ Heel Rise: _____ Swing Ext: _____ Stance Ext: _____ Stair/Ramp: _____

Obstacle Assist: On Off Sitting Mode: On Off Stance Lock: On Off

Activities: Cycling Golfing Sprinting Elliptical Rehabilitation Custom: _____

Notes: _____